

Dr. Han CMC Arthroplasty (Resection Trapeziectomy with LRTI)

GLOVES: 6.5 Micro Indicator/6.5 Micro

Anesthesia: General

PREP: Hibiclens

POSITION: SUPINE, Patient on Stretcher, Rolling Hand Table, Tourniquet, Extremity Drape, Sitting Stools

MEDICATIONS: 0.9% NAACL IRRIGATION 500 ML
20 mL, 0.25 Marcaine or 0.5% Marcaine

INSTRUMENTS: Hand Kit (Curved Littler or Small Curved Tenotomy Scissors, Adson Forceps), Drill, 2 and 3 mm Drill Bit Set, K-Wire Driver, **22 (or 24) Gauge Dental Wire, Freer Elevator, Cottle Elevator, Baby Holman Retractors, Heiss Retractor, Small Weitlaner Retractor, Mcglamry elevators (In Podiatry Set), Osteotomes, Ronguers (Wide Mouth), HEWSON SUTURE Retriever (If Available), A type of Small Cork Screw (ex. Arthrex AR-8919-01S)**

EQUIPMENT: MINI C-ARM (OrthoScan if available), Drill, K Wire Driver, Oscillating Sagittal Saw with a small short blade

SUPPLIES: HALF SHEET
Extremity Drape
Ray-Tech
LIGHT HANDLE COVER
15 BLADE (3)
Beaver Blade #0064 if not available #0069
10CC SYRINGE (1)
25 GAUGE HYPO (1)
BIPOLAR CORD/Forcep and Bipolar Foot Paddle
TOURNIQUET 18"
4" ESMARK (OPEN, STERILE)

SUTURES: 3-0 or 2-0 FiberLoop, 3-0 Ethibond x 2, and Straight Keith Needle
4-0 Vicryl RB-1, 3-0 Monocryl P-3, and 4-0 NYLON P-C 3

DRESSINGS: Steri Strips, Large Xeroform, WET 4X4, DRY 4X4, 4" Cast Padding x 2, 4" plaster x 12 sheet, 4" ACE WRAP, any kind of tape, SLING (IF HAD BLOCK),

OTHERS: Additional Equipment for Extra-Procedures **If Asked during Booking**
Arthrex Mini Tight Rope 1.1mm with C-Clamp Guide
Arthrex Corkscrew Micro and Nano Anchor (1st Preference for MP Arthroplasty) or
Biomet Juggerknot Soft Anchor 1.0mm Mini x 2 or
DePuy Synthes Micro QuickAnchor Suture Anchor (3-0 Needle) x 2 (Least Preferred)

LAST UPDATED: 1/30/2019